

Szczotki z nylonowym włosiem, jednorazowego użytku
uchwyt - skręcony drut ze stali nierdzewnej

NR KAT.	ŚREDNICA [mm]	PRZYBLIŻONA DŁUGOŚĆ CZĘŚCI CZYSZCZĄCEJ [mm]	PRZYBLIŻONA DŁUGOŚĆ CAŁKOWITA -BEZ PĘTLI [mm]	OPAK. [szt.]
IL.BR06039	1,0	6,35	152	2
IL.BR08039	1,0	12,7	203	2
IL.BR18039	1,0	6,35	457	2
IL.BR06062	1,57	12,7	152	2
IL.BR18062	1,57	12,7	457	2
IL.BR08070	1,77	12,7	203	2
IL.BR18070	1,77	12,7	457	2
IL.BR24070	1,77	12,7	610	2
IL.BR08079	2,0	19,0	203	2
IL.BR12079	2,0	19,0	305	2
IL.BR24080	2,0	25,4	610	2
IL.BR08087	2,2	12,7	203	2
IL.BR18087	2,2	12,7	457	2
IL.BR24098	2,49	12,7	610	2
IL.BR08105	2,67	25,4	203	2
IL.BR18105	2,67	25,4	457	2
IL.BR08118	3,0	38,1	203	2
IL.BR12118	3,0	69,8	305	2
IL.BR08131	3,3	19,0	203	2
IL.BR12158	4,0	19,0	305	2
IL.BR18158	4,0	19,0	457	2
IL.BR12197	5,0	76,2	305	2
IL.BR16197	5,0	76,2	406	2
IL.BR18197	5,0	76,2	457	2
IL.BR24197	5,0	76,2	610	2
IL.BR12236	6,0	50,8	305	2
IL.BR16236	6,0	50,8	406	2
IL.BR24236	6,0	50,8	610	2

Producent- Steris/Key Surgical

